

King, Valerie A.

From: Steve Hambek [Steve.Hambek@mdsps.com]
Sent: Thursday, April 29, 2004 9:35 AM
To: Kimberly Prchal
Subject: Re: Creatinine "Batches"

Sensitivity: Confidential

Kim,

The samples for each study (or group if applicable) are run on the same day.

Each day of operation is considered to be a batch/run. At the beginning of the day, standards and controls are run to verify performance of the assay. Controls are analyzed at intervals during the run to assure that the performance has not degraded.

The instrument is a "continuous feed" instrument so there is no "stop and start" that would constitute different runs.

Steven Hambek
Clinical Laboratory Manager, Lincoln
MDS Pharma Services